

 **AMD LASERS®**



**The Most Awarded  
Dental Lasers in the Industry**

World-Class Products | [amdlasers.com](http://amdlasers.com)

© 2023 AMD Lasers, Inc. All Rights Reserved.

# PICASSO

**THE FIRST DIODE LASER  
SPECIFICALLY DESIGNED  
FOR DENTAL HYGIENISTS**

*clarío*™

The Picasso Clario is the first diode laser specifically designed for dental hygienists and it sports a symbolic purple color. The streamlined navigation and pre-set parameters based on diagnosis facilitate simple, controlled, and appropriate use of the laser for hygiene treatments. Portable and affordable, the Picasso Clario laser puts the power and art of laser dentistry in the hands of the hygienist.

## FEATURES

- Easy to Use Presets
- Wireless Foot Control
- Built-In Treatment Timers for Periodontal Treatment and Aphthous Ulcer Treatment
- Multiple Languages Menu

## APPLICABLE PROCEDURES

- Hemostasis
- Coagulation
- Sicular Debridement
- Curettage
- And More...



## TECHNICAL SPECIFICATIONS

WEIGHT	DIMENSIONS	OUTPUT VOLTAGE	OUTPUT CURRENT	OUTPUT POWER RANGE	LASER CERTIFICATION
2.0 lbs (1 kg)	9.5" x 6" x 6.2"	+12.0 V DC	5.0 AMPS	60 Watts MAX	IV (4)
MEDIUM	WAVELENGTH	POWER MODES	DELIVERY FIBER	OUTPUT RANGE	
GaAlAs	810 ± 10nm	Continuous	400µm	0.5–1.2 Watts ±20% output	



VISIBLE AND INVISIBLE LASER RADIATION. AVOID EYE OR SKIN EXPOSURE TO DIRECT OR SCATTERED RADIATION.  
 WAVELENGTH: 810NM      WAVELENGTH: 850NM  
 CLASS 4 LASER PRODUCT      CLASS 2 LASER PRODUCT  
 OUTPUT: 1.2 WATTS      OUTPUT: 50mW  
 RAYONNEMENT LASER VISIBLE ET INVISIBLE. EXPOSITION DANGEREUSE DE L'ŒIL OU DE LA PEAU AU RAYONNEMENT DIRECT OU DIFFUS.  
 LA LONGUEUR D'ONDE: 810NM      APPAREIL À LASER DE CLASSE 4  
 ÉMISSION MAXIMALE: 1.2 WATTS      LA LONGUEUR D'ONDE: 850NM      APPAREIL À LASER DE CLASSE 2  
 ÉMISSION MAXIMALE: 50mW      IEC 60825-1:2014 - SAFETY OF LASER PRODUCTS

# PICASSO LITE+<sup>®</sup>

## THE MOST AWARDED DENTAL LASER IN THE INDUSTRY

The Picasso Lite+ Diode Laser was designed for versatility, usability, precision, and affordability. Offering a wide range of necessary capabilities for soft tissue surgical procedures and periodontal treatments, the Picasso Lite+ soft tissue laser is ideal for both dentists and hygienists alike. The Picasso Lite+ is equipped with a multi-tip handpiece, wireless foot control, and customizable modes & presets.

### FEATURES

- Easy to Use Presets
- Wireless Foot Control
- Built-In Treatment Timers for Periodontal Treatment and Aphthous Ulcer Treatment
- Multiple Languages Menu
- Continuous/ Pulsed Modes

### APPLICABLE PROCEDURES

- Frenectomy
- Fibroma Removal
- Troughing
- Hemostasis
- Periimplantitis
- Sicular Debridement
- Lip and Tongue Tie
- Tissue Tag
- Coagulation
- Curettage
- Herpetic Lesions
- And More...



### TECHNICAL SPECIFICATIONS

WEIGHT	DIMENSIONS	OUTPUT VOLTAGE	OUTPUT CURRENT	OUTPUT POWER RANGE	LASER CERTIFICATION
2.0 lbs (1 kg)	9.5" x 6" x 6.2"	+12.0 V DC	5.0 AMPS	60 Watts MAX	IV (4)
MEDIUM	WAVELENGTH	POWER MODES	DELIVERY FIBER	OUTPUT RANGE	
GaAlAs	810 ± 10nm	Continuous, Pulse Repeat	200µm - 400µm	0.5-3.0 Watts ±20% output	



VISIBLE AND INVISIBLE LASER RADIATION. AVOID EYE OR SKIN EXPOSURE TO DIRECT OR SCATTERED RADIATION.  
 WAVELENGTH: 810NM WAVELENGTH: 650NM  
 OUTPUT: 30mW OUTPUT: 30mW  
 CLASS 4 LASER PRODUCT CLASS 2 LASER PRODUCT  
 RAYONNEMENT LASER VISIBLE ET INVISIBLE. EXPOSITION DANGEREUSE DE L'ŒIL, OU DE LA PEAU AU RAYONNEMENT DIRECT OU DIFFUS.  
 LA LONGUEUR D'ONDE: 810NM APPAREIL À LASER DE CLASSE 4  
 ÉMISSION MAXIMALE: 3 WATTS APPAREIL À LASER DE CLASSE 2  
 ÉMISSION MAXIMALE: 30mW  
 IEC 60825-1:2014 - SAFETY OF LASER PRODUCTS

REV. 1

# PICASSO+

## THE COMPLETE DENTAL LASER SOLUTION

As AMD's top tier laser, the Picasso+ Diode Laser offers a complete solution for dental professionals looking for power, control, and flexibility. The high-power 7.0 wattage of the Picasso+ allows for laser-assisted whitening and a full range of technology such as a multi-tip handpiece, wireless foot control, customizable pulsed or continuous modes, and easy-to-use presets for a full range of soft-tissue and periodontal procedures.

### FEATURES

- Easy to Use Presets
- Wireless Foot Control
- Built-In Treatment Timers for Periodontal Treatment and Aphthous Ulcer Treatment
- Multiple Languages Menu
- Continuous/ Pulsed Modes

### APPLICABLE PROCEDURES

- Frenectomy
- Fibroma Removal
- Troughing
- Sicular Debridement
- Hemostasis
- Lip and Tongue Tie
- Perimplantitis
- Operculectomy
- Coagulation
- Teeth Whitening\*
- Curettage
- Tissue Tag
- Herpetic Lesions
- Implant Recovery
- And More...



### TECHNICAL SPECIFICATIONS

WEIGHT	DIMENSIONS	OUTPUT VOLTAGE	OUTPUT CURRENT	OUTPUT POWER RANGE	LASER CERTIFICATION
3.0 lbs (1.3 kg)	9.5" x 6" x 6.2"	+12.0 V DC	5.0 AMPS	60 Watts MAX	IV (4)
MEDIUM	WAVELENGTH	POWER MODES	DELIVERY FIBER	OUTPUT RANGE	WHITENING MODE
GaAlAs	810 ± 10nm	Continuous, Pulse Repeat	200µm - 400µm	0.5-7.0 Watts ±20% output	Preset to 7.0 Watts 30 second timer



VISIBLE AND INVISIBLE LASER RADIATION. AVOID EYE OR SKIN EXPOSURE TO DIRECT OR SCATTERED RADIATION.  
WAVELENGTH: 810NM      APPAREIL À LASER  
OUTPUT: 7 WATTS      DE CLASSE 4  
CLASS 4 LASER PRODUCT      CLASS 4 LASER PRODUCT  
RAYONNEMENT LASER VISIBLE ET INVISIBLE. EXPOSITION DANGEREUSE DE L'ŒIL, DE LA PEAU AU RAYONNEMENT DIRECT OU DIFFUS.  
LA LONGUEUR D'ONDE: 810NM      APPAREIL À LASER  
ÉMISSION MAXIMALE: 7 WATTS      DE CLASSE 4  
LA LONGUEUR D'ONDE: 810NM      APPAREIL À LASER  
ÉMISSION MAXIMALE: 50W      DE CLASSE 2  
\*IC 60825-2014 - SAFETY OF LASER PRODUCTS


Our goal is to ensure that every dental professional and practice has access to state-of-the-art advancements in medical & dental technology, as well as educational resources and training in order to elevate the standards of patient care through technology and innovation.

Along with our award winning line of diode lasers, **AMDLASERS.COM** also has many other product innovations available.

Including the Monet<sup>®</sup>, the world's first handheld laser curing light, Sheer White!<sup>®</sup> teeth whitening system, and many other great products.

GET THE MOST OUT OF YOUR PICASSO WITH OUR  
**SUBSCRIBE & SAVE**

PROGRAM FOR **20% OFF**  
REPLACEMENT TIPS AT  
**AMDLASERS.COM**





## **IMPLEMENTING LASERS IS EASY!**

Soft tissue lasers are more affordable than you might think, and implementing them into your practice is very easy. When you purchase a laser from AMD you will get access to free laser training, live product and setup demos, access to our laser expert, and many other resources to make implementation and daily use of your laser a simple process.

## **BETTER HEALING IS BETTER RETENTION**

Diode lasers help reduce tissue damage, bleeding, tooth loss, and receding gums. The Picasso soft tissue lasers streamline the elimination of periodontal disease and minimize the need for local anesthetic, stitches, and periodontal packs.

Using soft-tissue diode lasers for surgical procedures in your dental practice increases revenue and retention by decreasing time in the dentist chair and expediting patient healing.

AMD's Diode Lasers has the right laser for your needs. They can perform a variety of procedures, from Laser hygiene to periodontal and other soft-tissue laser procedures, allowing you to expand the services you offer and stop outsourcing periodontal laser procedures.



**SCAN NOW**  
OPT-IN FOR PROMOTIONS,  
FREE CE, AND MORE.



© 2023 AMD Lasers, Inc. All Rights Reserved.

**866.999.2635**  
[www.amdlasers.com](http://www.amdlasers.com)