

## QUARTERLY PROMOTIONS

# **SMDRT***Rewards* Make it Count with Your Next Order!


SM>RTRewards

### Enroll now at Shop.BrasselerUSA.com/ SmartRewards

to unlock exclusive, members-only savings when you purchase the quality dental essentials you trust.

> APRIL-JUNE 2023

# INTRODUCING THE BRASSELER USA Preparation-To-Polish Restorative System





## You know when something is just right

- PATIENT COMFORT Clipped corners provide a gentler fit
- DURABILITY Durable enough to handle all your operatory activities
- ► IMAGE CLARITY Technology designed to produce crystal clear images
- PLUG & PLAY Hands-off integration with most practice management platforms

#### **GEM™** Sensor Kit

- 1 GEM Digital X-Ray Sensor
- 100 GEM Sensor Barriers
- 1 Positioning System

Retail: \$6,995.00



(Sensors shown at actual size)

Product		GEM10	GEM15	GEM20			
Order No.		5028383U0 5028384U0		5028385U <b>0</b>			
Sensor Technology			CMOS				
External	L	36.36	38.83	41.80			
Dimensions	W	24.53	29.63	30.48			
(mm)	Н	12.52 13.01		13.02			
Active Imaging Are	a (mm)	30.01 x 20.01 32.99 x 25.82		35.92 x 25.82			
Active Image Area	(px)	1539 x 1026	1692 x 1324	1842 x 1324			
Image Resolution			20+ visible lp/mm				
Image Pixel Size		19.5 x 19.5 µm					
Maximum Gray Lev	vels	16,384 (14 bit)					
Cable Length			1.9m				

BRASSELER





#### **BLU Microlaser Starter Kit**

- 1 BLU microlaser
- 2 Lithium-ion batteries
- 1 Wireless foot pedal
- 2 AA batteries for foot pedal
- 1 Base charger
- 1 Power supply
- 10 Disposable tips (5 each: 5mm, 7mm)
- 3 Pairs of safety glasses
- 1 Pack of Initiating paper

Retail: \$6,495.00 Special: \$5,525.00



# **DIAGNOSTIC INSTRUMENTS**

Buy 24 Instruments, Get 45% Off!\* Buy 12 Instruments, Get 35% Off!\* Buy 6 Instruments, Get 25% Off!\*



# **RESTORATIVE MATERIALS**

#### DualForce<sup>™</sup> Matrix System Complete Kit 209022-B

Retail: \$571.99

#### Special: \$428.99

- 1 Premolar Ring
- 1 Molar Ring
- 25 Ultra-Wrap Matrix Bands, 4.5mm
- 25 Ultra-Wrap Matrix Bands, 5.5mm
- 25 Ultra-Wrap Matrix Bands, 6.5mm
- 1 Band Punch Pliers
- 1 Ring Forceps

- 20 #1 Active-Wedges 20 #2 Active-Wedges
- 20 #3 Active-Wedges
- 20 #4 Active-Wedges
- 20 #5 Active-Wedges



#### **G5™ All-Purpose Desensitizer**

- · Eliminates post-operative sensitivity in posterior composite restorations
- · Works by coagulating plasma proteins within dentinal tubule fluid, which helps eliminate movement of the fluid-the root cause of dentinal sensitivity
- · Used after Max Etch and before MPa Max adhesive

**5mL Bottle** 226105-B Retail: \$88.99 ea.

**Buy 2 5mL Bottles** of Desensitizer, Get 1 5mL Bottle at No Additional Cost!

**10mL Bottle** 226110-B Retail: \$127.99 ea.

Buy 2 10mL Bottles of Desensitizer, Get 1 10mL Bottle at No Additional Cost!

BRASSELE





Accur8<sup>™</sup> Dual-Arch Impression Posterior Trays 50-Pack 047501-B Retail: \$124.99

Buy 2 Packs, Get 1 at No Additional Cost!



# **DIRECT RESTORATION**

Brasseler USA's recommended diamonds, carbides, and polishers for each step of a Direct Restoration.

## COMPOSITE

Amalgam Removal Diamond	Ā	Re	ialgam moval rbides				Composit Removal Diamonds				Decay Removal Carbide	
					Sabrecut®					ET <sup>®</sup> Series		
Us No.	-	Us	No.	-	1558SC		Us No.		-	DET6F	Us No.	4
Order No.	6813	Ord	er No.	H34	H031RS		Order No.		6830	134F	Order No.	H1
Head Size (1/10 mm)	014	Hea	ad Size (1/10 mm)	010	012		Head Size (1)	/ 10 mm)	014	014	Head Size (1/10 mm)	014
Cutting Length (mm)	1.5	Cut	ting Length (mm)	2.0	4.0		Cutting Leng	th (mm)	2.7	6.0	Cutting Length (mm)	-
Grit	$\bigcirc$	Blac	de		-		Grit				Blade	-
Shank Type	FG	Sha	ank Type	FG	FG		Shank Type		FG	FG	Shank Type	RA
Diamonds			Micro	Carbi	uco	H	<u> </u>	C-Series	M	<u> </u>		
Us No.	330D	557D	-	Us No.		4C	35C	245C	330C	557C		
Order No.	D232	D835	838M	Order N	lo.	H1	H2	H245	H7	H31		
Head Size (1/10 mm)	008	010	007	Head S	ize (1/10 mm)	014	010	008	008	010		
Cutting Length (mm)	2.0	4.2	2.5	Cutting	Length (mm)	-	0.9	2.8	2.0	3.9		
Grit	$\bigcirc$	$\bigcirc$	$\bigcirc$	Blade		-	-	-	-	-		
Shank Type	FG	FG	FG	Shank	Туре	FG	FG	FG	FG	FG		
Finishing Diamonds						Finishi Carbide	•	Safe				

Diamonds		<u> </u>	<u> </u>	<u> </u>	<u> </u>
	ET® S	ieries			
Us No.	DET6F	DET9F	-	DOS1F	-
Order No.	134F	135F	8888	8379	8392
Head Size (1/10 mm)	014	014	012	023	016
Cutting Length (mm)	6.0	9.0	8.0	4.2	5.0
Grit					
Shank Type	FG	FG	FG	FG	FG

Buy 20 Packs, Get 40% Off!\* Buy 10 Packs, Get 30% Off!\* Buy 5 Packs, Get 20% Off!\*

#### A A A A ET<sup>®</sup> Series Safe-End Us No. ET6 ET9 7901 – 7408/0S1 Order No. H134 H135 H246 H274 H379 Head Size (1/10 mm) 014 014 009 016 023 Cutting Length (mm) 6.0 9.0 3.6 3.5 4.2 Blade Shank Type FG FG FG FG FG

#### **FINISHING CARBIDE BLADE CONFIGURATION**

12 Blade

#### **DIAMOND GRIT SIZE**

- Coarse
- Medium
- Fine

# **DIRECT RESTORATION**

				0.11			Get 30% Off!*
VisionFlex™	order No. Wid WS25A		gth (mm)	Grit 🔴 🔵 Asst. (EF, F	Single/Dou (, M) Sing		Buy 5 Packs,
Gateway™	DS25A		147	Asst. (EF, F			Get 20% Off!*
			(		Buy 20	Packs, Get 40	% Off!
eparating Strip					•	Packs, Get 30	
					Buy 5 P	acks, Get 20%	o Off!
Order No. Thi 9816	ckness (mm) 0.05	Width (mm) L 6.0	ength (mm) 150				
0010	0.00	0.0	100				
Iniversal	- Action	atio	- Alastan	dib			
composite	CON	III	EFIS	III		$\bigcirc$	
olishers		Ĭ					
			gle-Patient Sterile		ЕР™	Versaflex™	Composipro™
rder No.	SFLBCO-1	SBFLCP-CO	SFLBF-2	SBFLCP-F	EP1SM	VFD1LG	D209001
ead Size (1/10 mm) hape	100 Feather Lite™	055 Feather Lite™ Cup	100 Feather Lite™	055 Feather Lite <sup>™</sup> Cup	80 Pop-On Disc	120 Snap-On Disc	060 Brush
rit	Medium	Medium	Fine	Fine	C, M, F, XF	SC, C, M, F	Medium/Fine
hank Type	RA	RA	RA	RA	EPM21 Mandrel	310.21 Mandrel	RA
	Nelse	_	marca	_			
	Selle.		- AND	-			
	-	Diacomp® M	ulti-Use			Packs, Get 40	
			FLDC5	DC3.21	•	Packs, Get 30	
rder No.	FLDC5M	DC3M.21			Buy 5 P	acks, Get 20%	o U <b>⊞</b> *
<b>rder No.</b> ead Size (1/10 mm)	<b>FLDC5M</b> 140	060	140	060			
			140 1.6	060 7.5			
ead Size (1/10 mm) utting Length (mm) hape	140	060 7.5 Cup	1.6 Feather Lite™				
ead Size $(1/_{10} \text{ mm})$ utting Length (mm)	140 1.6	060 7.5	1.6	7.5			

Forza<sup>™</sup> LS-BR Latch BRASSELER ZALS Retail: \$519.00





Forza<sup>™</sup> Pro 1:1 Latch

Retail: \$1,479.00

Special: \$1,257.15

BRASSELER

Buy 2 at Special Price, Get 1 at No Additional Cost!

Mix and match. Must be of equal or lesser value.

# **INDIRECT RESTORATION**

Brasseler USA's recommended diamonds, carbides, and polishers for each step of Crown & Bridge Restorations.

Buy 20 Packs, Get 40% Off!\* Buy 10 Packs, Get 30% Off!\* Buy 5 Packs, Get 20% Off!\*

## CROWN REMOVAL

Porcelain, Zirconia & Lithium Disilicat Removal Diamond	· · · · · · · · · · · · · · · · · · ·	PFM, Met Substruct Cutting C	ture		Depth Cutting Diamonds				Ē	Ē	
	DuraCut™							Winter Restora	tive Conce	pts	
Order No.	6856DC	Order No.		H34L	Order No.		WDLCONV	WDLEXT	WDL1.0	WDL1.5	WDL2.0
Head Size (1/10 mm)	016	Head Size (1)	/ 10 mm)	012	Head Size (1/10 n	nm)	030	034	016	017	017
Cutting Length (mm)	8.0	Cutting Leng		3.5 I	Depth Cut (mm)		0.5 / 0.7 / 0.9	0.8 / 1.0 / 1.2	1.0	1.5	2.0
Grit		Blade			Grit		$\bigcirc$	$\bigcirc$	$\bigcirc$	$\bigcirc$	$\bigcirc$
Shank Type	FG	Shank Type		RA	Shank Type		FG	FG	FG	FG	FG
Breaking Contact Diamonds			Occlu Reduc Diame		H			Margin Refinement Diamonds		<b>—</b>	
Order No.	5858	5850	Order I	No.	5379	5811			Winter	<sup>r</sup> Restorativ	e Concept
Head Size (1/10 mm)	014	012		Size (1/10 mm)	023	033		Order No.		WSMR1	-
Cutting Length (mm)	8.0	10.0		Length (mm)	4.2	4.2		Head Size (1/10 mm)		010	
Grit	0.0		Grit					Cutting Length (mm)		0.27	
Shank Type	FG	FG	Shank	Tyne	FG	FG		Tip Size (1/10 mm)	,	1.0	
onunit typo	1 G		onum	1)00	14	1 G		Grit			
								Shank Type		FG	
Preparation Diamonds	Multi-Radi	us	ļ	Å				Shank Type		FG	
		F		F s Multi-Padius	Å					FG	
Diamonds	F C	Winter Restorat				K		KS Series			KS6
Diamonds Order No.	WS6M1.0	Winter Restorat WS6M1.2	WS8M0.6	WS8M1.0	WS8M1.2		50 K	KS Series S1 KS2	KS	3	<b>KS6</b> 025
Diamonds Order No. Head Size (1/10 mm)	<b>WS6M1.0</b> 016	Winter Restorat WS6M1.2 018	<b>WS8M0.6</b> 014	<b>WS8M1.0</b> 018	<b>WS8M1.2</b> 020	0	<b>50 К</b> а 10 0	KS Series 51 KS2 12 014	<b>KS</b> : 01	<b>3</b> 6	025
Diamonds Order No. Head Size (1/10 mm) Cutting Length (mm)	<b>WS6M1.0</b> 016 6.0	Winter Restorat WS6M1.2 018 6.0	<b>W\$8M0.6</b> 014 8.0	<b>WS8M1.0</b> 018 8.0	<b>WS8M1.2</b> 020 8.0	0	<b>50 К</b> а 10 0	KS Series S1 KS2	KS	<b>3</b> 6	
Diamonds Order No. Head Size (1/10 mm) Cutting Length (mm) Tip Size (1/10 mm)	<b>WS6M1.0</b> 016 6.0 1.0	Winter Restorat           WS6M1.2           018           6.0           1.2	<b>WS8M0.6</b> 014 8.0 0.6	<b>WS8M1.0</b> 018 8.0 1.0	WS8M1.2 020 8.0 1.2	0:	<b>SO K</b> 3 10 0 .0 7	KS Series           51         KS2           12         014           .0         8.0	<b>KS</b> 01 8.0	<b>3</b> 6	025 9.0 —
	<b>WS6M1.0</b> 016 6.0	Winter Restorat WS6M1.2 018 6.0	<b>W\$8M0.6</b> 014 8.0	<b>WS8M1.0</b> 018 8.0	<b>WS8M1.2</b> 020 8.0	0:	<b>SO K</b> 10 0 .0 7 	KS Series           51         KS2           12         014           .0         8.0	<b>KS</b> 01 8.0	<b>3</b> 6	025 9.0
Diamonds Order No. Head Size (1/10 mm) Cutting Length (mm) Tip Size (1/10 mm) Grit	<b>WS6M1.0</b> 016 6.0 1.0	Winter Restorat           WS6M1.2           018           6.0           1.2	<b>WS8M0.6</b> 014 8.0 0.6 ●●●●	<b>WS8M1.0</b> 018 8.0 1.0 ● ● ●	<b>WS8M1.2</b> 020 8.0 1.2 ● ● ●	0:	<b>SO K</b> 10 0 .0 7 	KS Series 51 KS2 12 014 .0 8.0  G FG Modified Shoulder	KS3 014 8.0 –	3 6 6 ) 1 Modified Shoulder	025 9.0 -
Diamonds Order No. Head Size (1/10 mm) Cutting Length (mm) Tip Size (1/10 mm) Grit	<b>WS6M1.0</b> 016 6.0 1.0	Winter Restorat           WS6M1.2           018           6.0           1.2	<b>WS8M0.6</b> 014 8.0 0.6 ●●●●	<b>WS8M1.0</b> 018 8.0 1.0 ● ● ●	<b>WS8M1.2</b> 020 8.0 1.2 ● ● ●	0:	<b>SO K</b> 10 0 .0 7 	KS Series 51 KS2 12 014 .0 8.0  G FG Modified	KS3 014 8.0 –	3 6 ) i Modified	025 9.0 -
Diamonds Order No. Head Size (1/10 mm) Cutting Length (mm) Tip Size (1/10 mm) Grit Shank Type Preparation Diamonds	<b>WS6M1.0</b> 016 6.0 1.0	Winter Restorat           WS6M1.2           018           6.0           1.2	WS8M0.6 014 8.0 0.6 FG	<b>WS8M1.0</b> 018 8.0 1.0 ● ● ●	<b>WS8M1.2</b> 020 8.0 1.2 ● ● ●	0:	<b>SO K</b> 10 0 .0 7 	KS Series 51 KS2 12 014 .0 8.0  G FG Modified Shoulder	KS3 011 8.0 	3 6 6 ) 1 Modified Shoulder	025 9.0 -
Diamonds Order No. Head Size (1/10 mm) Cutting Length (mm) Tip Size (1/10 mm) Grit Shank Type Preparation	WS6M1.0 016 6.0 1.0 FG	Winter Restorat WS6M1.2 018 6.0 1.2 FG	WS8M0.6 014 8.0 0.6 € FG	WS8M1.0 018 8.0 1.0 ♥♥♥ FG	WS8M1.2 020 8.0 1.2 ♥♥♥ FG	0:	<b>SO K</b> 3 10 0 .0 7 - G F	KS Series 51 KS2 12 014 .0 8.0 - G FG Modified Shoulder	<b>KS</b> 011 8.0 FG	3 6 0 Modified Shoulder	025 9.0 -
Diamonds Order No. Head Size (1/10 mm) Cutting Length (mm) Tip Size (1/10 mm) Grit Shank Type Preparation Diamonds Order No.	WS6M1.0 016 6.0 1.0 FG FG 856	Winter Restorat           WS6M1.2           018           6.0           1.2           FG           FG	WS8M0.6 014 8.0 0.6 ► FG	WS8M1.0 018 8.0 1.0 ♥♥♥ FG 856	WS8M1.2 020 8.0 1.2 ♥♥♥ FG 856	0:	SO K3 10 0 .0 7  G F 856	KS Series 51 KS2 12 014 .0 8.0  G FG Modified Shoulder Modified Shoulder 847KR	KS3 014 8.0 − FG	3 6 0 Modified Shoulder H 47KR	025 9.0 -
Diamonds Order No. Head Size (1/10 mm) Cutting Length (mm) Tip Size (1/10 mm) Grit Shank Type Preparation Diamonds Order No. Head Size (1/10 mm)	WS6M1.0 016 6.0 1.0 FG FG 856 012	Winter Restorat           WS6M1.2           018           6.0           1.2           FG           FG	WS8M0.6 014 8.0 0.6 ► FG	WS8M1.0         018         8.0         1.0         ●         FG	WS8M1.2         020         8.0         1.2         ●      <	0:	SO K3 10 0 .0 7 − G F 856 021	KS Series 51 KS2 12 014 .0 8.0  G FG Modified Shoulder Shoulder Modified Shoulder Shoulder Modified Shoulder	KS3 010 8.0 FG	3 6 0 Modified Shoulder 47KR 018	025 9.0 -

# **SPECIALTY RESTORATION**

## LAMINATE VENEER

Depth Cutting Diamonds		
US No.	LVS1	LVS2
Order No.	834	834
Head Size (1/10 mm)	021	016
Cutting Length (mm)	6.0	6.0
Depth Cut (mm)	0.5	0.3
Grit	$\bigcirc$	$\bigcirc$
Shank Type	FG	FG

Preparation Diamonds		
US No.	LVS3	LVS4
Order No.	6844	6844
Head Size (1/10 mm)	016	014
Cutting Length (mm)	10.0/1.5	10.0/1.5
Depth Cut (mm)	-	-
Grit	Coarse/Fine	Coarse/Fine
Shank Type	FG	FG


#### Buy 20 Packs, Get 40% Off!\* Buy 10 Packs, Get 30% Off!\* Buy 5 Packs, Get 20% Off!\*

## INLAY/ONLAY



## **PERI-IMPLANTITIS**

Developed with Dr. Dennis Tarnow and Dr. Paul Rosen for treatment of peri-implantitis



#### **FINISHING CARBIDE BLADE CONFIGURATION** DuraCut<sup>™</sup> (DC) 12 Blade ○ 30 Blade

**DIAMOND GRIT SIZE** 



\*Mix and match

## HIGH SPEED HANDPIECES




## **INDIRECT RESTORATION**

## INTRA-ORAL ADJUSTING & POLISHING OF ALL CERAMIC MATERIALS

Adjusting Stone	
	LD Grinder™
Order No.	LD11CZR
Head Size (1/10 mm)	030
Cutting Length (mm)	7.5
Qty/Pack	5
Grit	Coarse
Shank Type	RA

#### LD Grinder™ Point

Retail: \$98.75/Pk Special: \$74.06/Pk



Buy 20 Packs, Get 40% Off!\* Buy 10 Packs, Get 30% Off!\* Buy 5 Packs, Get 20% Off!\*

Dialite® Intra-Oral Polishers				
Order No.	W16	W17	W18	FLW14
Head Size (1/10 mm)	040	060	100	140
Cutting Length (mm)	10.0	7.5	2.5	1.6
Qty/Pack	5	5	5	3
21 - RA - Right-Angle L	atch (22 mm)			
Porcelain				
Coarse	W16DCHP21	W17DCHP21	W18DCHP21	FLW14C.RA
Medium	W16DMHP21	W17DMHP21	W18DMHP21	FLW14M.RA
O Fine	W16DFHP.21	W17DFHP21	W18DFHP.21	FLW14F.RA
Zirconia				
Medium	W16MZR.RA	W17MZR.RA	W18MZR.RA	FLW14MZR.RA
Fine	W16FZR.RA	W17FZR.RA	W18FZR.RA	FLW14FZR.RA
Lithium Disilicate				
Medium	W16MLD.RA	W17MLD.RA	W18MLD.RA	FLW14MLD.RA
Fine	W16FLD.RA	W17FLD.RA	W18FLD.RA	FLW14FLD.RA

Recommended Speed for Points, Cups & Knife Edge: 5,000-7,000 rpm

Recommended Speed for Feather Lite: 4,000-6,000 rpm; Final Polish: 2,000-3,000 rpm

#### Intra-Oral Dialite<sup>®</sup> Points, Cups, & Knife-Edge Polishers

Retail: \$86.50/Pk

#### Intra-Oral Dialite<sup>®</sup> Feather Lite<sup>™</sup> Polishers Retail: \$68.10/Pk

Buy 15 Packs, Get 40% Off!\* Buy 10 Packs, Get 30% Off!\* Buy 5 Packs, Get 20% Off!\*

Sterile Universal Ceramic Polishers			CER	Ű
Order No.	SBPT	SBCP	SFLB	SBFLCP
Head Size (1/10 mm)	040	060	100	055
Qty/Pack	25	25	25	25
21 - RA - Right-Angle La	tch (22 mm)			
🔵 Medium	SBPTCER-1	SBCPCER-1	SFLBCER-1	SBFLCP-CE
Fine	SBPTF-2	SBCPF-2	SFLBF-2	SBFLCP-F

Recommended Speed: 4,000-6,000 rpm; Final Polish: 2,000-3,000 rpm

Buy 15 Packs, Get 40% Off!\* Buy 10 Packs, Get 30% Off!\* Buy 5 Packs, Get 20% Off!\* **Sterile BrioShine™ Points & Cups** Retail: \$88.45/Pk

**Sterile BrioShine™ Feather Lite™** Retail: \$115.65/Pk

**Sterile BrioShine™ Feather Lite™ Cups** Retail: \$95.10/Pk



### EXTRA-ORAL ADJUSTING & POLISHING OF ALL CERAMIC MATERIALS

Adjusting Stones			9	<b>LD Grinder™ Tap</b> Retail: \$33.80 e	a. Retail	<b>'inder'</b> : \$47.
Order No.	LD13	LD14	LD15	Special: \$25.35	ea. Speci	al: \$3
Head Size (1/10 mm)	040	050	150			
Cutting Length (mm)	10.0	13.0	3.0	LD Grinder <sup>™</sup> Wh		
Grit	Coarse, Medium	Coarse, Medium	Coarse, Medium	Retail: \$60.85 e	a.	
Shape	Flat-End Taper	Cylinder	Wheel	Special: \$46.74	ea.	
Shank Type	HP	HP	HP			
Dialite® Extra-Oral Polishers Order No.	H2	R17	R20	L26	FL26	
Head Size $(1/_{10} \text{ mm})$	040	170	200	260	260	
Cutting Length (mm)	13.0	2.5	3.0	2.0	1.6	
Qty/Pack	1	1	1	1	1	
11 - HP - Straight Han	apiece (44.5 mm)					
Porcelain						
Coarse	H2DCHP	R17DCHP	-	L26DCHP.11	FL26CHP	
Medium	H2DMHP	R17DMHP	-	L26DMHP.11	FL26MHP	
Fine	H2DFHP	R17DFHP	-	L26DFHP.11	FL26FHP	
Zirconia						
Medium	H2MZR.HP	R17MZR.HP	R20MZR.HP	L26MZR.HP	FL26MZR.HP	
Fine Fine	H2FZR.HP	R17FZR.HP	R20FZR.HP	L26FZR.HP	FL26FZR.HP	
Lithium Disilicate						
Medium     Fine	H2MLD.HP	R17MLD.HP	R20MLD.HP	L26MLD.HP	FL26MLD.HP	
	H2FLD.HP	R17FLD.HP	R20FLD.HP	L26FLD.HP	FL26FLD.HP	

Recommended Speed for Points, Wheel & Knife Edge: 8,000-10,000 rpm

Recommended Speed for Feather Lite: 4,000-10,000 rpm; Final Polish: 2,000-5,000 rpm

**Extra-Oral Dialite® Points** Retail: \$32.40 ea.

**Extra-Oral Dialite® Small Wheel Polishers** Retail: \$58.15 ea.

**Extra-Oral Dialite® Knife-Edge Polishers** Retail: \$82.70 ea.

**Extra-Oral Dialite®** Large Wheel & **Feather Lite™ Polishers** Retail: \$64.80 ea.

Buy 15 Packs, Get 40% Off!\* Buy 10 Packs, Get 30% Off!\* Buy 5 Packs, Get 20% Off!\*

STRAIGHT HANDPIECES



Forza<sup>™</sup> F1S Retail: \$799.00 Special: \$719.10

Buy 2 at Special Price, Get 1 at No Additional Cost! Mix and match. Must be of equal or lesser value.

Air



Forza<sup>™</sup> LS Straight Nosecone Retail: \$479.00


Buy 3, Get 1 at No Additional Cost! Mix and match. Must be of equal or lesser value.



\*Mix and match


# **ESR<sup>™</sup> CM RECIPROCATING FILES**

ESR CM utilizes the latest NiTi technology to provide clinicians with superior safety, efficiency and performance. ESR CM files are also designed to preserve more coronal tooth structure compared to WaveOne® Gold.



## ESR<sup>™</sup> CM TECHNIQUE

- 1. After straight line access, determine working length (WL) with a size #10 K-file.
- 2. Use the ESR CM #25 primary/gauging file to estimate the canal size by engaging in three light strokes.
  - a. If the file reaches the apex easily the canal is "Basic" (larger) and you should finish with the Large ESR CM #40 or Extra-Large ESR CM #50 file.
  - b. If the ESR CM #25 does not reach the apex after two series of three light strokes the canal is "Advanced" (tighter) and you need to use a #15 K-file or the ESR CM Scout file to reestablish a safe glidepath to WL. Next, finish with an ESR CM #30 (which actually has a smaller max diameter than the #25). For extremely tight canals consider finishing with the Small ESR CM #20.
- 3. Obturate with BC Sealer<sup>™</sup> and the matching BC Points<sup>™</sup> for the last file used to reach WL.




ESR CM files are provided in sterile blister packs for optimal convenience and safety. If sterility is compromised prior to treatment, the files can be sterilized without any loss of functionality or performance.

STERILE R	File/Tip ID	Tip Size	21 mm	25 mm	31 mm	ESR™ CM Files (3/Pack)
	$\triangle$ White	15C*	5029334U0	5029335U0	5029336U0	Buy 50 Packs, Get 40% Off!
BRASSELER	$\triangle$ White	15	5028076U0	5028084U0	5028092U0	Buy 25 Packs, Get 25% Off!
One Brasseler Boulevard Savannah, Georgia 31419 ESR Instructions for Use	🛆 Yellow	20	5029126U0	5029127U0	5029128U0	Buy 10 Packs, Get 15% Off!
	🔺 Red	25	5028078U0	5028086U0	5028094U0	
	A Blue	30	5028079U0	5028087U0	5028095U0	
	Black	40	5028081U0	5028089U0	5028097U0	*15C scout file uses a different austenite metallurgy, making it more robust in tight/
	🛆 Yellow	50	5028082U0	5028090U0	5028098U0	curved canals.

# NEW! XP-3D SHAPER+


Brasseler USA's next generation adaptive file designed to streamline your procedures and give you greater control and safety. With a new, improved Shaper Tip, this file allows for more efficient navigation, making it easier to guide around curves and keeping it centered in the canal.



The XP-3D Shaper+ features Adaptive Core Technology, which enables the central core of the file to move freely and adapt to the natural morphology of the canal. This helps with debris removal, making the process more efficient without occluding dental tubules. Additionally, the turbulence generated by the XP-3D Shaper+ enhances the penetration of irrigants into the dentinal tubules, improving overall disinfection of the canal.

TRADITIONAL XP-CORE FILES ADA Debris Accumulation/Compaction Free-Floa

XP-3D SHAPER+ ADAPTIVE CORE Free-Floating with Turbulance



Using the XP-3D Shaper+ can help reduce debris accumulation and compaction, while also reducing the stress applied to the canal wall. This innovative file transitions from ISO #15 to ISO #30, allowing for shorter procedure times and greater agility.

Experience the new design of the XP-3D Shaper+, which offers greater predictability, simplicity, and control. Discover the benefits of the next generation adaptive file, and take your procedures to the next level.

XP-3D <sup>™</sup> Shaper+ (3/Pack) Retail: \$85.50/Pk	File/Tip
Special: \$59.85/Pk	🔺 Blue

'Tip ID	Tip Size	21 mm (Stopper)	25 mm (Stopper)	31 mm (Stopper)
Blue	30	5028597U0	5028598U0	5028599U0



# TRICON® ALL-IN-ONE IRRIGATION SOLUTION

Unlike traditional irrigants or other advanced 2:1 solutions, Triton works differently by avoiding the use of EDTA and CHX altogether. The non-NaOCI components in Triton proactively dissolve the dentinal debris allowing for a lower concentration of NaOCI to be exposed to organic debris without as much buffering. Synergistic and simultaneous dissolution of organic and inorganic debris permits the clinician to use lower volumes of irrigation solution and ensure maximum clinical efficiency.

## TRITON VS. TRADITIONAL IRRIGATION

- Multi-functional single irrigation solution reduces chair time, procedural steps, and overall irrigation costs
- ✓ Simultaneous organic and inorganic debris removal
- Patent-pending composition allows for a lower concentration of NaOCI to be used while providing improved efficacy
- ✓ More effective at smear layer removal vs. NaOCI with EDTA
- ✓ Dissolves organic tissue up to 1.7x faster than traditional NaOCI

## COST-EFFECTIVE

EndoSequence

Not only will Triton reduce chair time by eliminating the need for multiple irrigation solutions and sterile water rinses, but it also reduces the number of single-use disposable syringes and tips you'll require. In some cases, Triton can cut your total irrigation cost per procedure in half.



Retail: \$199.95 Special: \$149.97

## BC SEALER<sup>™</sup> IS TURNING 15!

15 years of clinical success, and there's so much to celebrate!

#### EndoSequence<sup>™</sup> BC Sealer<sup>™</sup>

- 1 2g Syringe
- 15 BC Tips™
- 1 Safety Data Sheet (SDS)

Retail: \$200.00

Special Birthday Promotion: \$128.88 PLUS Take 35% Off all Other Bioceramic Products Use Promo Code: BCBDAY

# TotalVac<sup>™</sup> IRRIGATION SYSTEM

TotalVac addresses the negative and positive pressure suction and irrigation needs of clinicians performing root canal therapy procedures. The design allows for single-handed positive pressure irrigation with simultaneous suction. The TotalVac kit can be configured to perform both positive and negative pressure irrigation using the included tips, which vary in length and diameter.

TotalVac is a highly modular system which allows versatility and access to a multitude of needles of different lengths and diameters, based on your specific needs.

#### TotalVac™ Kit

- 25 High Volume Evacuation and Surgical Suction Adaptor
- 25 TotalVac Tubing
- 25 Plastic Suction Tip
- 25 25 Gauge Short Delivery Tip
- 25 30 Gauge Side Vented Long Delivery Tip

Retail: \$200.00 Special: \$163.00



#### TotalVac™ Tubing (20/Pack)

Retail: \$100.00/Pk Special: \$85.00/Pk

25GA Short Irrigation Tip (20/Pack)

Retail: \$18.00/Pk Special: \$15.30/Pk



## High Volume Evacuation and Surgical Suction Adaptor (20/Pack)

Retail: \$37.50/Pk Special: \$31.88/Pk

Triton & TotalVac Kit Bundle \$312.00 (Retail Value: \$399.95)





Visit **BrasselerUSADental.com/Discover-TotalVac** to read an article and watch an instructional video

presented by Dr. Allen Ali Nasseh, DDS, MMSc





6.0

Introducing Brasseler's Forza Pro series. This line of premium electric attachments includes a 1:5 High Speed, 1:1 Latch and new 1:1 FG, ideal for margin refinement! The Forza Pro series features a scratch-resistant, glass coating, improving both aesthetics and durability. The heat-blocking technology of the ceramic ball bearing design helps to minimize the risk of patient burns. Forza Pro also offers the longest warranty in the industry at 3 years!



- ► Compatible with all E-Type Electric Motors
- Scratch-Resistant Glass Coating
- ► 5-Port Water Spray
- Heat Blocker Technology
- Small Head Size
- 3-Year Warranty

#### Forza™ Pro 1:5

Fiber Optic
Head size: 14.5mm x 9.5mm
3-Year Warranty

Retail: \$1,679.00 Special: \$1,511.10

#### Forza<sup>™</sup> Pro 1:1 Latch

Fiber Optic
Head size: 12.4mm x 9.5mm
3-Year Warranty

Retail: \$1,449.00 Special: \$1,304.10


#### Forza™ Pro 1:1 FG

Fiber Optic
Head size: 13.8mm x 9.5mm
3-Year Warranty

Retail: \$1,449.00 Special: \$1,304.10

# FORZA<sup>™</sup> LS LOW SPEED SYSTEM

The Forza LS low speed handpiece series from Brasseler USA features a convenient Quick Coupler design which converts all Forza LS air motors into quick disconnect devices. The series offers durable, single-piece construction attachments and powerful air motors for all major low speed brands: KaVo®, Midwest®, NSK®, Star®, and W&H®.



KaVo® is a registered trademark of Kaltenbach & Voigt GmbH, Germany. Midwest® is a registered trademark of DENTSPLY International. NSK® is a registered trademark of Nakanishi, Inc. Star® is a registered trademark of DentalEZ, Inc. W&H® is a registered trademark of W&H Dentalwerk Bürmoos GmbH, Austria.



# FORZA<sup>™</sup> ELXT ELECTRIC MOTOR


BRASSELER

OK

FORZA

DRZA ELXT

1:5 % I



- Reduced Risk of Patient Burns with the advanced Safety Check Function
- Low Cost, Hassle-free installation with the unique split-module space saving design
- Lightweight and compact motor provides optimal balance and comfort

\*Data on file.


#### Forza™ ELXT

- · LED Fiber Optics
- •Torque: 4.2 Ncm
- $\cdot$  Speed Range: 100-200,000 rpm
- · Motor Weight: 59g
- · 3-Year Warranty

Retail: \$3,349.00



# **FORZA<sup>™</sup> SHORT SERIES ELECTRIC ATTACHMENTS**





#### **BRASSELER U.S.A. DENTAL, LLC**

One Brasseler Blvd. Savannah, GA 31419

Toll Free: 800.841.4522 Phone: 912.925.8525 Fax: 888.610.1937

Shop.BrasselerUSA.com



©2023 Brasseler USA. All rights reserved. Printed in the U.S.A.



B-5720-US-72M-04.23-RRD





## TRUST BRASSELER USA FOR FAST, ECONOMICAL REPAIR OF YOUR HANDPIECES

We Service All Major Handpiece Brands including: A-dec°/W&H°, Bien-Air°, KaVo°, Midwest°, NSK°, Star°

HANDPIECE

QUALITYGREAT PRICE

► FAST TURNAROUND

Not sure if we can repair your handpiece? Contact your sales representative or contact Brasseler USA directly! 1.800.841.4522 (ext. 4675)

A-dec<sup>®</sup> is a registered trademark of A-dec, Inc., USA. Bien-Air<sup>®</sup> is a registered trademark of Bien-Air Medical Technologies, Switzerland. Ka\0<sup>®</sup> is a registered trademark of Kaltenbach & Voigt GmbH, Germany. Midwest<sup>®</sup> is a registered trademark of DENTSPLY International. NSK<sup>®</sup> is a registered trademark of Nakanishi, Inc. Star<sup>®</sup> is a registered trademark of DentalEZ, Inc. W&H<sup>®</sup> is a registered trademark of W&H Dentalwerk Bürmoos GmbH, Austria.

Invoice or statement prices may reflect or be subjected to a bundled discount or rebate pursuant to purchase offer, promotion, or discount program. You must fully and accurately report to Medicare, Medicaid, Tricare and/or any other federal or State program, upon request by such program, the discounted price(s) or net price(s) for each invoiced item, after giving effect to any applicable discounts or rebates, which price(s) may differ from the extended prices set forth on your invoice. Accordingly, you should tetain your invoice and all relevant information for your records. It is your responsibility to review any agreements or other documents, including offers or promotions, applicable to the invoiced products/prices to determine if your purchase(s) are subject to a bundled discount or rebate. Any such discounts must be calculated pursuant to the terms of the applicable purchase offer, promotion, or discount program. Participation in a promotional discount program is only permissible in accordance with discount program rules. By participation in such program, you agree that, to your knowledge, your practice complies with the discount program requirements.